


**Preparation of Samples for Illumina Next Generation Sequencing**

Example:



Design NGS-primers by adding the respective index 1-12 to the 5'-end of your selection primer sequence. Please attach the same index at 5'-end of forward- as well as of reverse primer. It may be necessary to replace 5'-nucleotides from either or both of your original primers in order to keep annealing temperatures of both primers in the same range.

Indexed primers:

<b>Forward-Primer (5'→3')</b>	
FP-Index 1	<u>ATCACG</u> -FP-Sequence
FP-Index 2	<u>CGATGT</u> FP-Sequence
FP-Index 3	<u>TTAGGC</u> FP-Sequence
FP-Index 4	<u>TGACCA</u> FP-Sequence
FP-Index 5	<u>ACAGTG</u> FP-Sequence
FP-Index 6	<u>GCCAAT</u> FP-Sequence
FP-Index 7	<u>CAGATC</u> FP-Sequence
FP-Index 8	<u>ACTTGA</u> FP-Sequence
FP-Index 9	<u>GATCAG</u> FP-Sequence
FP-Index 10	<u>TAGCTT</u> FP-Sequence
FP-Index 11	<u>GGCTAC</u> FP-Sequence
FP-Index 12	<u>CTTGTA</u> FP-Sequence
FP-Index 13	<u>AGTCAA</u> FP-Sequence

Reverse-Primer (5'→3')		
RP-Index 1	<u>ATCACG</u>	RP-Sequence
RP-Index 2	<u>CGATGT</u>	RP-Sequence
RP-Index 3	<u>TTAGGC</u>	RP-Sequence
RP-Index 4	<u>TGACCA</u>	RP-Sequence
RP-Index 5	<u>ACAGTG</u>	RP-Sequence
RP-Index 6	<u>GCCAAT</u>	RP-Sequence
RP-Index 7	<u>CAGATC</u>	RP-Sequence
RP-Index 8	<u>ACTTGA</u>	RP-Sequence
RP-Index 9	<u>GATCAG</u>	RP-Sequence
RP-Index 10	<u>TAGCTT</u>	RP-Sequence
RP-Index 11	<u>GGCTAC</u>	RP-Sequence
RP-Index 12	<u>CTTGTA</u>	RP-Sequence
RP-Index 13	<u>AGTCAA</u>	RP-Sequence

- PCR using Taq Core Kit from Qiagen:

- PCR-Mix (for 100 µl-PCR)\*

	1x
10x PCR-Puffer*:	10
5'-Primer-Index (100 µM):	1
3'-Primer-Index (100 µM):	1
dNTPs (10 mM)*:	2
Templat (approx: 1µM)	2
H <sub>2</sub> O:	83,5
Taq*:	0,5
Total (µl)	100

\* Please calculate the number of PCR reactions that are needed to end up with the necessary amount of DNA (depends on number of differently indexed amplicons derived from different experiments that have to be combined to give one sequencing probe. For one sequencing probe minimum of xµg of dsDNA is needed.

- PCR-program:

1.	94°C	60sec	(denaturation)
2.	xx °C	60sec	(annealing)
3.	72°C	60sec	(elongation)
4.	GOTO 1	rep 5-8	
5.	8°C	endless	

- Check on 2.5% agarose
- Purification with QiaQuick PCR Purification Kit or MiniElute PCR Purification Kit (Qiagen)
- Concentration and purity measurement of each indexed probe (A260/A280)
- Pooling of all indexed probes (sum: minimum 3 µg) (optional: please inquire AptaIT if pooling should be done in your case)
- Shipment of the probes according to the guidelines of the NGS-provider (refrigerated/chilled in EB or H<sub>2</sub>O).